

Laboratory Press TP 20, TP 40, TP 60

Manually operated hydraulic laboratory press

TP20 TP40 TP60

This manually operated hydraulic press allows all the compacting operations common in the laboratory to be carried out.

A special press tool available as an accessory is used to produce pellets.

The hydraulic pump is operated by a hand lever. The direction of movement of the plunger is reversed by switching over the valve to permit simple access to the dies installed.

The threaded spindle serves as a counter-bearing to absorb the compacting pressure and to reduce the no-load stroke.

We reserve the right to make technical changes.



print

