

Pelletsing Press TP 40/2d, TP 60/2d

Manually operated hydraulic pelletizing press

TP40/2d TP60/2d

By means of this manually operated hydraulic pelletizing press it is possible to produce very easy tablets with different diameters (depending on the spectrometer).

For this purpose a special press tool with the corresponding diameter is inserted into the press and can be filled upon moving back the upper cross beam. The press tool is delivered as special accessory.

The hydraulic pump is operated by the handle. The direction of motion of the press piston can be reserved by changing over the valve. The threaded spindle is intended for counter piece of the pressure absorption and for the reduction of the idling stroke.

We reserve the right to make technical changes.



print

